



**Product Specification**  
**Faeces tube, with spoon, push cap, (LxØ): 75 x 23.5 mm, transparent**

**Page 1**



**Product description**

Order number 80.620  
Product description Faeces tube, with spoon, transparent, push cap, cap: brown, (LxØ) without cap: 75 x 23.5 mm, 50 piece(s)/bag

**Product characteristics**

Closure type push cap  
Design spoon

**Size**

Diameter 23.5 mm  
Length of product 87 mm  
Length excluding cap 75 mm

**Material & colours**

Colour of product transparent  
Closure material Low Density Polyethylene (LD-PE)  
Colour of cap brown  
Tube material Polystyrene (PS)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Faeces tube, with spoon, push cap, (LxØ): 75 x 23.5 mm, transparent**

**Purity & certification**

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Batched          | yes                           |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 500           |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 50            |
| Piece (s) / outer case        | 500           |
| Piece(s) / pallet             | 18000         |
| Depth of case                 | 398 mm        |
| Width of case                 | 398 mm        |
| Height of case                | 300 mm        |
| Case volume                   | 0.0475 cbm    |
| Weight of product             | 0.0091 kg     |
| Weight of case                | 5.29 kg       |
| EAN of inner box              | 4038917082952 |
| EAN case                      | 4038917001939 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

