



Product Specification

Faeces tube, with spoon, screw cap, (LxØ): 107 x 25 mm, transparent, sterile

Page 1



Product description

| | |
|---------------------|--|
| Order number | 80.622 |
| Product description | Faeces tube, with spoon, transparent, screw cap, Screw cap with stool spoon holder, cap: brown, (LxØ) without cap: 107 x 25 mm, sterile, 50 piece(s)/bag |

Product characteristics

| | |
|--------------|-----------------------------------|
| Cap | Screw cap with stool spoon holder |
| Closure type | screw cap |
| Design | spoon |
| Base shape | flat-conical base |

Size

| | |
|----------------------|--------|
| Diameter | 25 mm |
| Length of product | 111 mm |
| Length excluding cap | 107 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-11-15

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Faeces tube, with spoon, screw cap, (LxØ): 107 x 25 mm, transparent, sterile

Page 2

Material & colours

| | |
|-------------------|-----------------------------------|
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | brown |
| Stem material | Polystyrene (PS) |
| Tube material | Polypropylene (PP) |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 250 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 250 |
| Piece(s) / pallet | 10000 |
| Depth of case | 598 mm |
| Width of case | 398 mm |
| Height of case | 165 mm |
| Case volume | 0.0393 cbm |
| Weight of product | 0.0129 kg |
| Weight of case | 4.29 kg |
| EAN of inner box | 4038917082976 |
| EAN case | 4066916010805 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-11-15

This document was prepared by EDP support and is valid without signature.



SARSTEDT