



Product Specification

Faeces tube, with special faeces spoon (1 ml), screw cap, (LxØ): 101 x 16.5 mm, transparent

Page 1



Product description

| | |
|---------------------|---|
| Order number | 80.623.022 |
| Product description | Faeces tube, with special faeces spoon (1 ml), With enclosed stirring rod, transparent, screw cap, cap: brown, (LxØ) without cap: 101 x 16.5 mm, with paper label, 250 piece(s)/bag |

Product characteristics

| | |
|-----------------------|-----------------------------|
| Label/ Print | with paper label |
| Colour of print/label | black |
| Closure type | screw cap |
| Design | special faeces spoon (1 ml) |

Size

| | |
|----------------------|---------|
| Diameter | 16.5 mm |
| Length of product | 107 mm |
| Length excluding cap | 101 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-11-15

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Faeces tube, with special faeces spoon (1 ml), screw cap, (LxØ): 101 x 16.5 mm, transparent

Page 2

Material & colours

| | |
|-------------------|-----------------------------------|
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | brown |
| Tube material | Polypropylene (PP) |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 250 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 24000 |
| Depth of case | 398 mm |
| Width of case | 398 mm |
| Height of case | 415 mm |
| Case volume | 0.0657 cbm |
| Weight of product | 0.0079 kg |
| Weight of case | 8.62 kg |
| EAN of inner box | 4038917115407 |
| EAN case | 4038917014601 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-11-15

This document was prepared by EDP support and is valid without signature.



SARSTEDT