



**Product Specification**  
**Faeces tube, with spoon, screw cap, (LxØ): 76 x 20 mm, transparent**



**Product description**

Order number	80.734
Product description	Faeces tube, with spoon, transparent, screw cap, Screw cap with stool spoon, assembled, cap: brown, (LxØ) without cap: 76 x 20 mm, 100 piece(s)/bag

**Product characteristics**

Cap	Screw cap with stool spoon
Closure type	screw cap
Design	spoon
Base shape	flat-conical base

**Size**

Diameter	20 mm
Length of product	77 mm
Length excluding cap	76 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Faeces tube, with spoon, screw cap, (LxØ): 76 x 20 mm, transparent**

**Page 2**

**Material & colours**

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown
Stem material	Polystyrene (PS)
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	21000
Depth of case	388 mm
Width of case	378 mm
Height of case	235 mm
Case volume	0.0345 cbm
Weight of product	0.007 kg
Weight of case	4.07 kg
EAN of inner box	4038917083812
EAN case	4038917002707

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**