



Product description

| | |
|---------------------|--|
| Order number | 93.853.134 |
| Product description | Block Rack D13, for 50 collection tubes, Ø opening: 13 mm, format: 5 x 10, blue, material: recycled PP |

Product characteristics

| | |
|------------------------|------------|
| Label/ Print | with print |
| Number of wells | 50 |
| Storage capacity | 50 |
| Format | 5 x 10 |
| Autoclavable (max. °C) | 121 °C |

Size

| | |
|---------------------|--------|
| Diameter of opening | 13 mm |
| Width of product | 96 mm |
| Height of product | 45 mm |
| Length of product | 180 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material recycled Polypropylene (recycled PP)
Colour of product blue

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest subpackaging -
Piece(s) / inner box 1
Piece (s) / outer case 16
Piece(s) / pallet 1440
Depth of case 398 mm
Width of case 368 mm
Height of case 120 mm
Case volume 0.0176 cbm
Weight of product 0.1687 kg
Weight of case 3.16 kg
EAN of inner box 4038917122559
EAN case 4038917074056

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

