



## Product Specification

Medicine cup, max. volume: 50 ml, (LxØ): 54 x 55 mm, graduated, PP, transparent

Page 1



### Product description

Order number	75.9924.522
Product description	Medicine cup, max. volume: 50 ml, (LxØ): 54 x 55 mm, Ø opening: 55.4 mm, transparent, graduated, material: PP

### Product characteristics

Graduation	yes
Application	Administration of medicines
Example of use	Oral administration of medications in liquid or tablet form
Volume of work	50 ml
Autoclavable (max. °C)	121 °C

### Size

Sample volume	50 ml
Max. volume	50 ml
Diameter	55 mm
Diameter of opening	55,4 mm
Height of product	54 mm
Length of product	54 mm
Length excluding cap	55.4 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

Medicine cup, max. volume: 50 ml, (LxØ): 54 x 55 mm, graduated, PP, transparent

Page 2

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

### Purity & certification

Product category	Medical device
CE certified	CE
Latex-free	yes
Pyrogen-free	no
Safety data sheet	no

### Packaging

Minimum order qty.	500
Piece(s) / inner box	500
Piece (s) / outer case	500
Piece(s) / pallet	80000
Depth of case	271 mm
Width of case	224 mm
Height of case	158 mm
Case volume	0.0099 cbm
Weight of case	3 kg
EAN of inner box	4038917352949
EAN case	4038917352932

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**