



### **Product description**

Order number	95.1706
Product description	Wall holder, material: stainless steel, for Formalin system, 1 piece(s)/case

### **Product characteristics**

Application	Tissue fixation
Suitable for	Formalin system

### **Size**

Width of product	101 mm
Length of product	200 mm

### **Material & colours**

Product material	stainless steel
------------------	-----------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      yes

### **Packaging**

Minimum order qty.                      1  
Type of smallest                              case  
subpackaging  
Piece(s) / inner box                      1  
Piece (s) / outer case                      1  
Piece(s) / pallet                              324  
Depth of case                              268 mm  
Width of case                                191 mm  
Height of case                               100 mm  
Case volume                                0.0051 cbm  
Weight of product                         0.269 kg  
Weight of case                              0.42 kg  
EAN of inner box                         4038917335461  
EAN case                                    4038917335454

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

