



### **Product description**

Order number	73.1055
Product description	RIA tube, suitable for gamma counters, material: PP, 0.5 ml, (LxØ): 22 x 8 mm, transparent, 1,000 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	push cap
Suitable for	gamma counters
Base shape	conical base

### **Size**

Max. volume	0.5 ml
Diameter	8 mm
Diameter of opening	29 mm
Height of product	22 mm
Length of product	22 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**RIA tube, suitable for gamma counters,**  
**transparent**

**Page 2**

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Batched no

**Packaging**

Minimum order qty. 10000  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 1000  
Piece (s) / outer case 10000  
Piece(s) / pallet 480000  
Depth of case 388 mm  
Width of case 378 mm  
Height of case 215 mm  
Case volume 0.0315 cbm  
Weight of product 0.0005 kg  
Weight of case 5.05 kg  
EAN of inner box 4038917085175  
EAN case 4038917003865

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-09-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**