



### **Product description**

Order number	73.680
Product description	Mini vial, cap enclosed, material: HD-PE, 6 ml, (LxØ): 56 x 15 mm, white, 1,000 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	screw cap
Base shape	flat base

### **Size**

Max. volume	6 ml
Diameter	15 mm
Diameter of opening	13 mm
Height of product	56 mm
Length of product	56 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	High Density Polyethylene (HD-PE)
Colour of product	white
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	white

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	48000
Depth of case	570 mm
Width of case	215 mm
Height of case	370 mm
Case volume	0.0453 cbm
Weight of product	0.0033 kg
Weight of case	7.4 kg
EAN of inner box	4038917083348
EAN case	4038917002288

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

