



Product description

Order number	62.610.018
Product description	Screw cap tube, working volume: 10 ml, (LxØ): 92 x 15.3 mm, material: PP, round base, transparent, screw cap, natural, cap assembled, with print, label/print: green, "CSF/Liquor", with graduations, sterile, 1 piece(s)/blister

Product characteristics

Label/ Print	with print
Colour of print/label	green
Graduation	yes
Closure type	screw cap
Base shape	round base
Volume of work	10 ml
Centrifugation max (RCF)	3000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Screw cap tube, 10 ml, (LxØ): 92 x 15.3 mm, PP,
with print**

Page 2

Size

Diameter 15.3 mm
Length of product 92 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap natural
Tube material Polypropylene (PP)

Purity & certification

Satisfies the requirement IATA
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest blister
subpackaging
Piece(s) / inner box 100
Piece (s) / outer case 500
Piece(s) / pallet 13500
Depth of box 351 mm
Width of box 151 mm
Height of box 171 mm
Depth of case 778 mm
Width of case 368 mm
Height of case 195 mm
Case volume 0.0558 cbm
Weight of product 0.0052 kg
Weight of case 4.23 kg
EAN of inner box 4038917233132
EAN case 4038917233101

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT