



Product description

Order number 93.1650.010
Product description Insertion plate, Version 1, material: PP, 1 piece(s)/case

Product characteristics

Suitable for Safety Tray (93.1650) & e.g. 1x Multi-Safe twin plus 500 (77.3884.005/.105), 1x rack with handle (93.851 - 93.851.1xx), 2x dispenser boxes (95.813.001), 3x multi-purpose cup (75.1355), 1x disinfection spray

Size

Width of product 272 mm
Length of product 420 mm

Material & colours

Product material Polypropylene (PP)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest
subpackaging case
Piece(s) / inner box 1
Piece (s) / outer case 1
Piece(s) / pallet 1440
Case volume 0.0009 cbm
Weight of product 0.152 kg
Weight of case 0.31 kg
EAN of inner box 4038917132695
EAN case 4038917076715

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

