



## Product description

Order number	77.3884.305
Product description	Disposal container, Multi-Safe twin plus REC, 500 ml, grey, with plastic label, label/print: yellow/red, biohazard labelling, material: recycled PP, application e.g.: disposal of pointed and sharp objects, opening: round, height: 245 mm, Ø of cap: 82 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 126 piece(s)/carton

## Product characteristics

Label/ Print	with plastic label
Colour of print/label	yellow/red
Example of use	disposal of pointed and sharp objects
Absorbent liner	yes
Shape of opening	round
Removal aid - turning aid	Removal and rotation aid

## Size

Max. volume	500 ml
Diameter	116 mm
Diameter of base	105 mm
Diameter of lid	82 mm
Diameter of opening	71 mm
Height of product	245 mm

## Material & colours

Product material	recycled Polypropylene (recycled PP)
Colour of product	grey
Closure material	Polypropylene (PP)
Colour of cap	red
Colour of lid	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Satisfies the requirement	TRBA 250, UN 3291, DIN EN ISO 23907
Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	126
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	126
Piece(s) / pallet	3024
Depth of case	599 mm
Width of case	266 mm
Height of case	451 mm
Case volume	0.0719 cbm
Weight of product	0.0753 kg
Weight of case	10.36 kg
EAN of inner box	4066916004163
EAN case	4066916004170

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

