



## Product description

Order number	77.3884.350
Product description	Disposal container, Multi-Safe twin plus REC, 5,000 ml, grey, with plastic label, label/print: yellow/red, biohazard labelling, material: recycled PP, application e.g.: disposal of pointed and sharp objects, opening: oval, height: 379 mm, lid (LxW): 179 x 176 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 34 piece(s)/carton

## Product characteristics

Label/ Print	with plastic label
Colour of print/label	yellow/red
Example of use	disposal of pointed and sharp objects
Absorbent liner	yes
Shape of opening	oval
Removal aid - turning aid	Removal and rotation aid

## Size

Max. volume	5000 ml
Diameter of base	112 mm
Height of product	379 mm

## Material & colours

Product material	recycled Polypropylene (recycled PP)
Colour of product	grey
Closure material	Polypropylene (PP)
Colour of cap	red
Colour of lid	red

## Purity & certification

Satisfies the requirement	TRBA 250, UN 3291, DIN EN ISO 23907
Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Disposal container, Multi-Safe twin plus REC,**  
**5,000 ml, biohazard labelling**

**Page 2**

**Packaging**

Minimum order qty.	34
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	34
Piece(s) / pallet	408
Depth of case	800 mm
Width of case	400 mm
Height of case	380 mm
Case volume	0.1216 cbm
Weight of product	0.342 kg
Weight of case	13.11 kg
EAN of inner box	4066916004446
EAN case	4066916004453

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-09-01

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**