



Product description

Order number	77.3890.203
Product description	Disposal container, Multi-Safe mini REC, 300 ml, grey, material: recycled PP, opening: round, height: 160 mm, lid (LxW): 81 x 45 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 135 piece(s)/carton

Product characteristics

Shape of opening	round
Removal aid - turning aid	Removal aid

Size

Max. volume	300 ml
Diameter of base	80 mm x 45 mm
Diameter of opening	27 mm
Height of product	160 mm

Material & colours

Product material	recycled Polypropylene (recycled PP)
Colour of product	grey
Colour of cap	red
Colour of lid	red

Purity & certification

Satisfies the requirement	TRBA 250, UN 3291, DIN EN ISO 23907
Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	135
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	135
Piece(s) / pallet	2700
Depth of case	600 mm
Width of case	400 mm
Height of case	310 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Disposal container, Multi-Safe mini REC, 300 ml

Page 2

Case volume	0.0744 cbm
Weight of product	0.0694 kg
Weight of case	10 kg
EAN of inner box	4066916004569
EAN case	4066916004576

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT