



### **Product description**

Order number	95.901
Product description	Post transport packaging, application e.g.: external packaging according to P 650/ ADR, high design, (LxWxH): 107 x 198 x 50 mm, for diagnostic specimens, holds 1 to 4 transport containers/bottles up to a diameter of 30 mm, Locking mechanism, 50 piece(s)/case
UN identification	yes

### **Product characteristics**

Label/ Print	with print
Colour of print/label	black
Design	high design
Example of use	external packaging according to P 650/ ADR
Suitable for	diagnostic specimens

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Width of product	198 mm
Height of product	50 mm
Length of product	107 mm

### **Material & colours**

Product material	cardboard
Colour of product	brown

### **Purity & certification**

Satisfies the requirement	P 650, ADR, IATA
Product category	no medical device   IVD
Batched	yes

### **Packaging**

Minimum order qty.	50
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	50
Piece(s) / pallet	2100
Depth of case	656 mm
Width of case	414 mm
Height of case	115 mm
Case volume	0.0312 cbm
Weight of product	0.0701 kg
Weight of case	4.04 kg
EAN of inner box	4038917084567
EAN case	4038917064583

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

