




Product Specification

Cold transport container for blood gas capillaries,
S-Monovette® up to 105 x Ø 18 mm, transparent,
length: 50 mm

Page 1



Product description

Order number	95.997
Product description	Cold transport container for blood gas capillaries, S-Monovette® up to 105 x Ø 18 mm, transparent, length: 50 mm, Cold shipping / warm shipping, 1 piece(s)/case

Product characteristics

Application	Cold shipping / warm shipping
Suitable for	blood gas capillaries, S-Monovette® up to 105 x Ø 18 mm
Absorbent liner	no

Size

Width of product	50 mm
Height of product	140 mm
Length of product	50 mm
Usable length	118 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Cold transport container for blood gas capillaries,
S-Monovette® up to 105 x Ø 18 mm, transparent,
length: 50 mm

Page 2

Material & colours

Product material	High Density Polyethylene (HD-PE)
Colour of product	transparent
Colour of cap	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	66
Piece(s) / pallet	2508
Depth of case	568 mm
Width of case	433 mm
Height of case	120 mm
Case volume	0.0295 cbm
Weight of product	0.0446 kg
Weight of case	4.59 kg
EAN of inner box	4038917084918
EAN case	4038917362382

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT