



Product description

Order number	95.1411
Product description	Mailing envelope, material: PE, blue, (LxW): 310 x 123 mm, 50 piece(s)/case

Product characteristics

Absorbent liner	no
-----------------	----

Size

Width of product	123 mm
Length of product	310 mm

Material & colours

Product material	Polyethylene (PE)
Colour of product	blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Satisfies the requirement	P 650
Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	1000
Piece(s) / pallet	40000
Depth of case	428 mm
Width of case	323 mm
Height of case	130 mm
Case volume	0.018 cbm
Weight of product	0.0081 kg
Weight of case	8.42 kg
EAN of inner box	4038917086547
EAN case	4038917005098

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

