



Product description

Order number	65.676.001
Product description	Screw cap, brown, suitable for mailing container 85 x 30 mm, 50 piece(s)/bag

Product characteristics

Application	Mailing of samples
Suitable for	mailing container 85 x 30 mm

Size

Length of product	15 mm
-------------------	-------

Material & colours

Product material	Polypropylene (PP)
Colour of product	brown

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 500
Type of smallest
subpackaging bag
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 39000
Depth of case 418 mm
Width of case 278 mm
Height of case 130 mm
Case volume 0.0151 cbm
Weight of product 0.0028 kg
Weight of case 1.71 kg
EAN of inner box 4038917115797
EAN case 4038917014885

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

