



### **Product description**

Order number	65.679
Product description	Screw cap, natural, suitable for Mailing container 126 x 30 mm, 50 piece(s)/bag

### **Product characteristics**

Application	Mailing of samples
Suitable for	Mailing container 126 x 30 mm

### **Size**

Length of product	20 mm
-------------------	-------

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	natural

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

**Packaging**

Minimum order qty.                      250  
Type of smallest  
subpackaging                              bag  
Piece(s) / inner box                      50  
Piece (s) / outer case                      250  
Piece(s) / pallet                            41250  
Depth of case                              338 mm  
Width of case                                228 mm  
Height of case                               120 mm  
Case volume                                0.0092 cbm  
Weight of product                         0.0035 kg  
Weight of case                              1.1 kg  
EAN of inner box                          4038917083331  
EAN case                                      4038917002271

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

