



### **Product description**

Order number	65.679.004
Product description	Screw cap, yellow, suitable for Mailing container 126 x 30 mm, 50 piece(s)/bag

### **Product characteristics**

Application	Mailing of samples
Suitable for	Mailing container 126 x 30 mm

### **Size**

Length of product	20 mm
-------------------	-------

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	yellow

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category no medical device | IVD  
Batched no

**Packaging**

Minimum order qty. 250  
Type of smallest subpackaging bag  
Piece(s) / inner box 50  
Piece (s) / outer case 250  
Piece(s) / pallet 41250  
Depth of case 338 mm  
Width of case 228 mm  
Height of case 120 mm  
Case volume 0.0092 cbm  
Weight of product 0.0035 kg  
Weight of case 1.1 kg  
EAN of inner box 4038917115834  
EAN case 4038917014939

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

