



## Product Specification

**Mailing container, transparent, construction: round, with absorbent liner, length: 85 mm, Ø opening: 30 mm, without cap**

**Page 1**



### Product description

Order number	78.892
Product description	Mailing container, round, transparent, construction: round, with absorbent liner, length: 85 mm, Ø opening: 30 mm, without cap, cap: natural, 50 piece(s)/bag

### Product characteristics

Application	Mailing of samples
Example of use	Secondary packaging according to P 650
Absorbent liner	yes
Shape of opening	round
Liquid-tight	yes

### Size

Diameter of opening	30 mm
Length of product	85 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Mailing container, transparent, construction:  
round, with absorbent liner, length: 85 mm, Ø  
opening: 30 mm, without cap**

**Page 2**

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Colour of cap	natural

### Purity & certification

Satisfies the requirement	P 650, ADR
Product category	no medical device   IVD
Batched	no

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	12000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0097 kg
Weight of case	5.69 kg
EAN of inner box	4038917084499
EAN case	4038917003308

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-06-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**