



Product Specification
**Urine container NFT, 3 l, with integrated
needle-free transfer unit, brown, with light
protection, graduated**

Page 1



Product description

Order number	77.575.900
Product description	Urine container NFT, Urine collection with hygienically closed, needle-free collection function for Urine Monovettes, max. working volume: 3 l, Ø opening: 83.5 mm, brown, with light protection, graduated, material: PE, with paper label, german/english, cap: white, with integrated needle-free transfer unit, assembled, 30 piece(s)/case

Product characteristics

Label/ Print	with paper label
Colour of print/label	white/black
Graduation	yes
Cap	membrane cap
Closure type	screw cap
Light protection	yes
Design	with integrated needle-free transfer unit
Application	Urine collection with hygienically closed, needle-free collection function for Urine Monovettes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**Urine container NFT, 3 l, with integrated
needle-free transfer unit, brown, with light
protection, graduated**

Page 2

Example of use	closed urine transfer into a Monovette® Urine
Volume of work	3.000 ml
Viewing window	yes
Connection	Luer (m)

Size

Sample volume	3.000.00 ml
Max. volume	3 l
Diameter of opening	83,5 mm
Width of product	130 mm
Length of product	125 mm
Length including cap	246 mm
Length excluding cap	243 mm

Material & colours

Product material	Polyethylene (PE)
Colour of product	brown
Closure material	Polypropylene (PP)
Colour of cap	white

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Safety data sheet	no
Batched	yes

Packaging

Minimum order qty.	30
Type of smallest subpackaging	case
Piece(s) / inner box	30
Piece (s) / outer case	30
Piece(s) / pallet	270
Depth of case	643 mm
Width of case	403 mm
Height of case	515 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**Urine container NFT, 3 l, with integrated
needle-free transfer unit, brown, with light
protection, graduated**

Page 3

Case volume	0.1335 cbm
Weight of product	0.1433 kg
Weight of case	5.1 kg
EAN of inner box	4038917947718
EAN case	4038917947725

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT