



### **Product description**

Order number	77.577
Product description	Wide neck bottle, Urine collection, max. working volume: 0.25 l, Ø opening: 35 mm, white, material: HD-PE, cap: black, 120 piece(s)/case

### **Product characteristics**

Graduation	no
Closure type	screw cap
Application	Urine collection
Viewing window	no

### **Size**

Sample volume	250 ml
Max. volume	0.25 l
Diameter of opening	35 mm
Width of product	60 mm
Height of product	90 mm
Length of product	60 mm
Length including cap	110 mm
Length excluding cap	98 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	High Density Polyethylene (HD-PE)
Colour of product	white
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	black

### **Purity & certification**

Product category	no medical device   IVD
Safety data sheet	no
Batched	yes

### **Packaging**

Minimum order qty.	120
Type of smallest subpackaging	case
Piece(s) / inner box	120
Piece (s) / outer case	120
Piece(s) / pallet	2880
Depth of case	458 mm
Width of case	288 mm
Height of case	415 mm
Case volume	0.0547 cbm
Weight of product	0.0371 kg
Weight of case	5.03 kg
EAN of inner box	4038917082624
EAN case	4038917001694

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

