



Product description

| | |
|---------------------|--|
| Order number | 75.1337.502 |
| Product description | Medicine cup, max. volume: 30 ml, (LxØ): 40 x 37 mm, Ø opening: 37 mm, red, graduated, material: PP, 75 piece(s)/bag |

Product characteristics

| | |
|----------------|-------|
| Graduation | yes |
| Volume of work | 30 ml |

Size

| | |
|----------------------|-------|
| Sample volume | 30 ml |
| Max. volume | 30 ml |
| Diameter | 37 mm |
| Diameter of opening | 37 mm |
| Height of product | 40 mm |
| Length of product | 40 mm |
| Length excluding cap | 40 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Medicine cup, max. volume: 30 ml, (LxØ): 40 x 37 mm, graduated, PP, red

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product red

Purity & certification

Product category Medical device
CE certified CE supplier
Batched yes

Packaging

Minimum order qty. 3750
Type of smallest subpackaging bag
Piece(s) / inner box 3750
Piece (s) / outer case 3750
Piece(s) / pallet 180000
Depth of case 190 mm
Width of case 390 mm
Height of case 390 mm
Case volume 0.0289 cbm
Weight of product 0.0015 kg
Weight of case 6 kg
EAN of inner box 4038917127738
EAN case 4042836002038

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT