



## Product Specification

Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, PP, with safety label, with integrated transfer unit, transparent

Page 1



### Product description

|                     |  |
|---------------------|--|
| Order number        | 75.562.400   |
| Product description | Container with screw cap, collection and storage of urine, max. working volume: 100 ml, (ØxH): 57 x 76 mm, Ø opening: 62 mm, transparent, graduated, material: PP, screw cap, cap: yellow, with safety label, with integrated transfer unit, cap assembled, sterile, 5 piece(s)/tube bag |

### Product characteristics

|                       |   |
|-----------------------|---|
| Label/ Print          | with safety label                               |
| Colour of print/label | white/red/black                                 |
| Graduation            | yes   |
| Closure type          | screw cap                                       |
| Design                | with integrated transfer unit                   |
| Application           | collection and storage of urine                 |
| Example of use        | closed urine transfer into a V-Monovette® Urine |
| Volume of work        | 20 ml - 100 ml                                  |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, PP, with safety label, with integrated transfer unit, transparent

Page 2

### Size

|                      |        |
|----------------------|--------|
| Sample volume        | 100 ml |
| Max. volume          | 100 ml |
| Diameter             | 57 mm  |
| Diameter of opening  | 62 mm  |
| Height of product    | 76 mm  |
| Length excluding cap | 72 mm  |

### Material & colours

|                   |                                   |
|-------------------|-----------------------------------|
| Product material  | Polypropylene (PP)                |
| Colour of product | transparent                       |
| Closure material  | High Density Polyethylene (HD-PE) |
| Colour of cap     | yellow                            |

### Purity & certification

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Purity standard  | sterile                       |
| Sterilisation    | Electron irradiation          |
| Batched          | yes                           |

### Packaging

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 200           |
| Type of smallest subpackaging | tube bag      |
| Piece(s) / inner box          | 200           |
| Piece (s) / outer case        | 200           |
| Piece(s) / pallet             | 5600          |
| Depth of case                 | 588 mm        |
| Width of case                 | 398 mm        |
| Height of case                | 255 mm        |
| Case volume                   | 0.0597 cbm    |
| Weight of product             | 0.0173 kg     |
| Weight of case                | 4.31 kg       |
| EAN of inner box              | 4038917236805 |
| EAN case                      | 4038917236799 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**