



Product description

Order number	75.563
Product description	Container with screw cap, collection and storage of urine, max. working volume: 100 ml, Ø: 57 mm, Ø opening: 57 mm, transparent, graduated, material: PP, screw cap, without cap, 50 piece(s)/bag

Product characteristics

Graduation	yes
Closure type	screw cap
Application	collection and storage of urine
Volume of work	100 ml

Size

Sample volume	100 ml
Max. volume	100 ml
Diameter	57 mm
Diameter of opening	57 mm
Length of product	72 mm
Length excluding cap	72 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 500
Type of smallest bag
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 15000
Depth of case 443 mm
Width of case 323 mm
Height of case 330 mm
Case volume 0.0472 cbm
Weight of product 0.0084 kg
Weight of case 5.01 kg
EAN of inner box 4038917082501
EAN case 4038917001588

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

