



### **Product description**

Order number	76.564
Product description	Screw cap, suitable for 100-ml container 75.563, material: HD-PE, natural, 50 piece(s)/bag

### **Product characteristics**

Closure type	screw cap
Example of use	cap for urine collection tube
Suitable for	100-ml container 75.563

### **Size**

Diameter	62 mm
Height of product	14 mm

### **Material & colours**

Product material	High Density Polyethylene (HD-PE)
Colour of product	natural

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	21000
Depth of case	388 mm
Width of case	378 mm
Height of case	250 mm
Case volume	0.0367 cbm
Weight of product	0.0068 kg
Weight of case	3.91 kg
EAN of inner box	4038917082518
EAN case	4038917001595

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

