



Product description

| | |
|---------------------|--|
| Order number | 60.548.012 |
| Product description | Screw cap tube, working volume: 50 ml, (LxØ): 114 x 28 mm, material: PP, conical base, transparent, HD-PE, without cap, with graduations, 250 piece(s)/bag |

Product characteristics

| | |
|--------------------------|--------------|
| Graduation | yes |
| Closure type | screw cap |
| Base shape | conical base |
| Volume of work | 50 ml |
| Centrifugation max (RCF) | 20000 x g |

Size

| | |
|-------------------|--------|
| Diameter | 28 mm |
| Length of product | 114 mm |

Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |

Purity & certification

| | |
|---------------------------|--------------------------------|
| Satisfies the requirement | ADR |
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Pyrogen-free | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|--------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 250 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 10000 |
| Depth of case | 598 mm |
| Width of case | 398 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP

Page 2

| | |
|-------------------|---------------|
| Height of case | 350 mm |
| Case volume | 0.0833 cbm |
| Weight of product | 0.0098 kg |
| Weight of case | 5.72 kg |
| EAN of inner box | 4038917111973 |
| EAN case | 4038917057868 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT