



### Product description

Order number	10.252.028
Product description	Urine Monovette® , nominal volume: 10 ml, without preparation, Urine tube, (LxØ): 102 x 15.3 mm, transparent, graduated, material: PP, Luer screw cap, cap yellow, with plastic label, sterile, 100 piece(s)/bag

### Product characteristics

Type of collection	Aspiration
Type of preparation	without preparation
Label/ Print	with plastic label
Colour of print/label	transparent/yellow
Graduation	yes
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Design	Luer lock
Application	Urine tube
Example of use	Urine collection, transport, and storage
Base shape	round base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine Monovette® , 10 ml, cap yellow, (LxØ): 102 x 15.3 mm, 100 piece(s)/bag**

**Page 2**

Volume of work	10 ml
Recommended centrifugation	400 x g
Centrifugation time	5 Min
Connection	Luer lock

**Size**

Sample volume	10 ml
Diameter	15.3 mm
Length including cap	117 mm
Length excluding cap	102 mm
Length including cap and plunger	127 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

**Purity & certification**

Satisfies the requirement	IATA, ADR
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	100
Piece(s) / pallet	7000
Case volume	0.0229 cbm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Urine Monovette® , 10 ml, cap yellow, (LxØ): 102 x 15.3 mm, 100 piece(s)/bag**

**Page 3**

Weight of product	0.0149 kg
Weight of case	1.54 kg
EAN of inner box	4038917104494
EAN case	4038917075671

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**