



Product Specification

**Urine Monovette® , Boric acid, 3.2 ml, cap green,
(LxØ): 75 x 13 mm, 64 piece(s)/bag**

Page 1

Product description

Order number	10.256.001
Product description	Urine Monovette® , Boric acid, nominal volume: 3.2 ml, Boric acid, Urine tube, (LxØ): 75 x 13 mm, transparent, fill mark, material: PP, Luer screw cap, cap green, incl. collection tip, with paper label, 64 piece(s)/bag

Product characteristics

Type of collection	Aspiration
Type of preparation	Boric acid
Label/ Print	with paper label
Colour of print/label	green
Graduation	fill mark
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Base shape	round base
Volume of work	3,2 ml
Recommended centrifugation	400 x g
Centrifugation time	5 Min
Connection	Luer (m)

Size

Sample volume	3.2 ml
Diameter	13 mm
Length including cap	90 mm
Length excluding cap	75 mm
Length including cap and plunger	109 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	green
Tube material	Polypropylene (PP)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-12-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Urine Monovette® , Boric acid, 3.2 ml, cap green,
(LxØ): 75 x 13 mm, 64 piece(s)/bag

Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Safety data sheet yes
Batched yes

Packaging

Minimum order qty. 512
Type of smallest bag
subpackaging
Piece(s) / inner box 64
Piece (s) / outer case 512
Piece(s) / pallet 36864
Depth of case 388 mm
Width of case 358 mm
Height of case 150 mm
Case volume 0.0208 cbm
Weight of product 0.0054 kg
Weight of case 3.1 kg
EAN of inner box 4038917982825
EAN case 4038917982832

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

