



## Product description

Order number	10.258.061
Product description	Urine Monovette® , Z, without preparation, Urine tube, (LxØ): 92 x 15 mm, transparent, graduated, material: PP, Luer screw cap, cap yellow, incl. collection tip, with paper label, suitable for Tempus600, 64 piece(s)/bag

## Product characteristics

Type of collection	Aspiration
Type of preparation	without preparation
Label/ Print	with paper label
Colour of print/label	yellow/white
Graduation	yes
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Suitable for	Tempus600
Base shape	round base
Volume of work	8,5 ml
Recommended centrifugation	400 x g
Centrifugation time	5 Min
Connection	Luer (m)

## Size

Diameter	15 mm
Length including cap	106 mm
Length excluding cap	92 mm
Length including cap and plunger	127 mm

## Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Tube material	Polypropylene (PP)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine Monovette® , Z, cap yellow, (LxØ): 92 x 15 mm, 64 piece(s)/bag**

**Page 2**

Piston material High Density Polyethylene (HD-PE)  
Piston plunger material Polystyrene (PS)

**Purity & certification**

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Safety data sheet no  
Batched yes

**Packaging**

Minimum order qty. 512  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 64  
Piece (s) / outer case 512  
Piece(s) / pallet 24576  
Case volume 0.0315 cbm  
Weight of product 0.0078 kg  
Weight of case 4.49 kg  
EAN of inner box 4038917992435  
EAN case 4038917992442

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-08-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**