



Product description

Order number	10.258.088
Product description	Urine Monovette® , nominal volume: 8.5 ml, without preparation, Urine tube, (LxØ): 92 x 15 mm, transparent, graduated, material: PP, Luer screw cap, cap yellow, incl. collection tip, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix, 64 piece(s)/bag

Product characteristics

Type of collection	Aspiration
Type of preparation	without preparation
Label/ Print	with plastic label pre-barcoded
Colour of print/label	white/yellow
Graduation	yes
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Base shape	round base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Urine Monovette® , 8.5 ml, cap yellow, (LxØ): 92 x 15 mm, 64 piece(s)/bag

Page 2

Volume of work	8,5 ml
Recommended centrifugation	400 x g
Centrifugation time	5 Min
Connection	Luer (m)

Size

Sample volume	8.5 ml
Diameter	15 mm
Length including cap	106 mm
Length excluding cap	92 mm
Length including cap and plunger	127 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitrodiagnostic device, CE
CE certified	CE
Safety data sheet	no
Batched	yes

Packaging

Minimum order qty.	20480
Type of smallest subpackaging	bag
Piece(s) / inner box	64
Piece (s) / outer case	512
Piece(s) / pallet	24576
Depth of case	388 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Urine Monovette® , 8.5 ml, cap yellow, (LxØ): 92 x 15 mm, 64 piece(s)/bag

Width of case	378 mm
Height of case	215 mm
Case volume	0.0315 cbm
Weight of product	0.0078 kg
Weight of case	4.52 kg
EAN of inner box	4038917951302
EAN case	4038917951319

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

