



Product description

Order number	11.1240.100
Product description	Urine transfer unit for V-Monovette® Urine, for Urine transfer, 50 piece(s)/minigrip bag

Product characteristics

Type of collection	Vacuum
Light protection	no
Design	With wings on the tube holder
Application	Urine transfer
Suitable for	V-Monovette® Urine
Connection	Vacuum tube

Size

Length of tube	10 cm
----------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Colour of product yellow

Purity & certification

Product category In vitrodiagnostic device, CE

CE certified CE supplier

Safety data sheet no

Batched yes

Packaging

Minimum order qty. 500

Type of smallest subpackaging minigrip bag

Piece(s) / inner box 50

Piece (s) / outer case 500

Piece(s) / pallet 30000

Depth of box 300 mm

Width of box 250 mm

Height of box 10 mm

Depth of case 380 mm

Width of case 350 mm

Height of case 170 mm

Case volume 0.0226 cbm

Weight of product 0.002 kg

Weight of case 1.8 kg

EAN of inner box 6947462401157

EAN case 6947462401423

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

