




Product Specification
Seraplas® valve filter, white, for separation of serum/plasma from the blood cells after centrifugation



Product description

Order number 53.677
Product description Seraplas® valve filter, V11, white, for separation of serum/plasma from the blood cells after centrifugation, 500 piece(s)/bag

Product characteristics

Application venous blood collection
Suitable for separation of serum/plasma from the blood cells after centrifugation

Size

Length of product 92 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Seraplas® valve filter, white, for separation of serum/plasma from the blood cells after centrifugation

Page 2

Material & colours

Product material Polyethylene (PE)
Colour of product white

Purity & certification

Product category no medical device | IVD
Safety data sheet no
Batched no

Packaging

Minimum order qty. 5000
Type of smallest subpackaging bag
Piece(s) / inner box 500
Piece (s) / outer case 5000
Piece(s) / pallet 120000
Depth of case 398 mm
Width of case 398 mm
Height of case 415 mm
Case volume 0.0657 cbm
Weight of product 0.0015 kg
Weight of case 8.5 kg
EAN of inner box 4038917083317
EAN case 4038917002257

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT