



Product description

Order number	55.1571.051
Product description	Adapter tube, (LxØ): 65 x 13 mm, material: PP, light blue, 250 piece(s)/bag

Product characteristics

Suitable for	S-Monovette® with Ø 8 mm, Microvette® and micro tubes
Base shape	round base

Size

Diameter	13 mm
Length of product	65 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	light blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1000
Type of smallest
subpackaging bag
Piece(s) / inner box 250
Piece (s) / outer case 1000
Piece(s) / pallet 72000
Depth of case 333 mm
Width of case 303 mm
Height of case 190 mm
Case volume 0.0192 cbm
Weight of product 0.0019 kg
Weight of case 2.23 kg
EAN of inner box 4038917283380
EAN case 4038917283373

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

