



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 65.647.020   |
| Product description | Archiving cap, with filter, light blue, suitable for S-Monovette® , tubes Ø 13-16 mm, 1,000 piece(s)/bag |

### **Product characteristics**

|              |                                 |
|--------------|---------------------------------|
| Suitable for | S-Monovette® , tubes Ø 13-16 mm |
| Filter       | yes                             |

### **Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | light blue         |
| Colour of cap     | light blue         |

### **Purity & certification**

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Batched          | yes                           |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Archiving cap, light blue, suitable for**  
**S-Monovette® , tubes Ø 13-16 mm**

**Page 2**

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 5000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 1000          |
| Piece (s) / outer case        | 5000          |
| Piece(s) / pallet             | 210000        |
| Depth of case                 | 388 mm        |
| Width of case                 | 378 mm        |
| Height of case                | 250 mm        |
| Case volume                   | 0.0367 cbm    |
| Weight of product             | 0.001 kg      |
| Weight of case                | 5.37 kg       |
| EAN of inner box              | 4038917250085 |
| EAN case                      | 4038917250078 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**