



Product Specification
Archiving cap, light blue, suitable for
S-Monovette® , tubes Ø 13-16 mm

Page 1



Product description

Order number	65.647.300
Product description	Archiving cap, without filter, light blue, suitable for S-Monovette® , tubes Ø 13-16 mm, 1,000 piece(s)/bag

Product characteristics

Suitable for	S-Monovette® , tubes Ø 13-16 mm
Filter	no

Material & colours

Product material	Polypropylene (PP)
Colour of product	light blue

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Archiving cap, light blue, suitable for
S-Monovette[®], tubes Ø 13-16 mm

Page 2

Packaging

Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	210000
Depth of case	378 mm
Width of case	250 mm
Height of case	388 mm
Case volume	0.0367 cbm
Weight of case	5.37 kg
EAN of inner box	4038917988315
EAN case	4038917988322

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT