



Product description

Order number	14.1112
Product description	Membrane adapter, connection: Luer, for S-Monovette® needle/Multifly® needle, 1 piece(s)/blister, 100 piece(s)/case, sterile, non-pyrogenic/endotoxin-free
PZN number	6796187

Product characteristics

Example of use	Adaptation of an S-Monovette® Needle or (Safety) Multifly® Needle to a Luer system, e.g. Blood-Gas Monovette®
Suitable for Connection	S-Monovette® needle/Multifly® needle Luer

Size

Length of product	16 mm
-------------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material High Density Polyethylene (HD-PE)
Colour of product translucent

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile, non-pyrogenic/endotoxin-free
Sterilisation Electron irradiation
Pyrogen-free yes
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest blister
subpackaging
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 80000
Depth of box 268 mm
Width of box 81 mm
Height of box 71 mm
Depth of case 373 mm
Width of case 178 mm
Height of case 287 mm
Case volume 0.0191 cbm
Weight of product 0.0005 kg
Weight of case 1.49 kg
EAN of inner box 4038917085359
EAN case 4038917004053

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

