



### Product description

Order number	51.1638.940
Product description	Safety-Multifly® BKF set 21G, needle (ØxL): 21G x 3/4", coding colour: green, tube length: 200 mm, with blood culture adapter set, set components: 1 Safety-Multifly® needle 21G with long tubing, DEHP-free (REF 85.1638.235) and 1 Blood Culture Adapter Universal (REF 14.1207), 1 piece(s)/blister, sterile, 1 piece(s)/case

### Product characteristics

Design	blood culture adapter set
Coding colour	green
Safety product	yes
Connection	blood culture bottle

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Safety-Multifly® BKF set 21G, 21G x 3/4", green,**  
**tube length: 200 mm, 1 piece(s)/blister**

**Page 2**

**Size**

Diameter of needle (external)	21 G / 0,8 mm
Diameter of tube	2.3 mm
Length of needle	19 mm
Length of tube	200 mm

**Purity & certification**

Product category	Medical device
Purity standard	sterile
Batched	yes

**Packaging**

Minimum order qty.	50000
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	50
Piece(s) / pallet	3200
Depth of case	398 mm
Width of case	298 mm
Height of case	215 mm
Case volume	0.0255 cbm
Weight of product	0.0082 kg
Weight of case	1.09 kg
EAN of inner box	4038917408738
EAN case	4038917408721

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-09-01

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**