



**Product Specification**  
**Safety-Multifly® needle, 23G x 3/4", blue, tube**  
**length: 240 mm, 1 piece(s)/blister**

**Page 1**



**Product description**

|                     |                                                                                                                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 85.1640.201                                                                                                                                                                                         |
| Product description | Safety-Multifly® needle, needle (ØxL): 23G x 3/4", coding colour: blue, tube length: 240 mm, with Luer (f), DEHP-free, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 120 piece(s)/case |

**Product characteristics**

|                |          |
|----------------|----------|
| Design         | Luer (f) |
| Coding colour  | blue     |
| Safety product | yes      |
| Connection     | Luer (f) |

**Size**

|                               |               |
|-------------------------------|---------------|
| Diameter of needle (external) | 23 G / 0,6 mm |
| Diameter of tube              | 2.3 mm        |
| Length of needle              | 19 mm         |
| Length of tube                | 240 mm        |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Safety-Multifly® needle, 23G x 3/4", blue, tube**  
**length: 240 mm, 1 piece(s)/blister**

**Page 2**

**Purity & certification**

|                  |                                       |
|------------------|---------------------------------------|
| Product category | Medical device                        |
| CE certified     | CE 0124 (DEKRA)                       |
| Purity standard  | sterile, non-pyrogenic/endotoxin-free |
| Sterilisation    | Gas sterilisation                     |
| Pyrogen-free     | yes                                   |
| Plasticiser      | DEHP-free                             |
| Batched          | yes                                   |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 480           |
| Type of smallest subpackaging | blister       |
| Piece(s) / inner box          | 120           |
| Piece (s) / outer case        | 480           |
| Piece(s) / pallet             | 20640         |
| Depth of box                  | 224 mm        |
| Width of box                  | 214 mm        |
| Height of box                 | 171 mm        |
| Depth of case                 | 442 mm        |
| Width of case                 | 358 mm        |
| Height of case                | 240 mm        |
| Case volume                   | 0.038 cbm     |
| Weight of product             | 0.003 kg      |
| Weight of case                | 3.32 kg       |
| EAN of inner box              | 4038917423946 |
| EAN case                      | 4038917423656 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**