




**Product Specification**  
Sample tube, EDTA K3E, 10 ml, cap red, (LxØ): 95 x 16.8 mm, with print

Page 1



### Product description

Order number	26.358
Product description	Sample tube, nominal volume: 10 ml, cap red, preparation: K3 EDTA, push cap, (LxØ): 95 x 16.8 mm, round base, with print, 100 piece(s)/bag

### Product characteristics

Type of preparation	K3 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with print
Colour of print/label	red
Graduation	yes
Cap	Skirted push cap
Closure type	push cap
Base shape	round base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Sample tube, EDTA K3E, 10 ml, cap red, (LxØ): 95  
x 16.8 mm, with print**

**Page 2**

**Size**

Sample volume	10 ml
Diameter	16.8 mm
Length of product	95 mm
Length excluding cap	95 mm

**Material & colours**

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	red
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	30000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.006 kg
Weight of case	6.63 kg
EAN of inner box	4038917081306
EAN case	4038917000444

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**