



**Product Specification**  
**Sample tube, Serum CAT, 7 ml, cap black, (LxØ):**  
**100 x 13 mm**

**Page 1**



**Product description**

Order number	29.373
Product description	Sample tube, nominal volume: 7 ml, cap black, preparation: Clotting activator, push cap, (LxØ): 100 x 13 mm, round base, 100 piece(s)/bag

**Product characteristics**

Type of preparation	Clotting activator
Colour code	EU
Graduation	no
Closure type	push cap
Base shape	round base
Recommended centrifugation	2000 x g
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-06-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Sample tube, Serum CAT, 7 ml, cap black, (LxØ):**  
**100 x 13 mm**

**Size**

Sample volume	7 ml
Diameter	13 mm
Length of product	100 mm
Length excluding cap	100 mm

**Material & colours**

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	black
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	42000
Depth of case	388 mm
Width of case	378 mm
Height of case	250 mm
Case volume	0.0367 cbm
Weight of product	0.0047 kg
Weight of case	5.22 kg
EAN of inner box	4038917081368
EAN case	4038917335966

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

