



**Product Specification**  
**Sample tube, Serum CAT, 4.5 ml, cap red, (LxØ):  
75 x 13 mm, with paper label**

**Page 1**



### Product description

Order number	32.329.100
Product description	Sample tube, nominal volume: 4.5 ml, cap red, preparation: Clotting activator, push cap, (LxØ): 75 x 13 mm, round base, with paper label, 100 piece(s)/bag

### Product characteristics

Type of preparation	Clotting activator
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	red
Graduation	fill mark
Closure type	push cap
Base shape	round base
Recommended centrifugation	2000 x g
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-11-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Sample tube, Serum CAT, 4.5 ml, cap red, (LxØ):**  
**75 x 13 mm, with paper label**

**Page 2**

**Size**

Sample volume	4.5 ml
Diameter	13 mm
Length of product	75 mm
Length excluding cap	75 mm

**Material & colours**

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	red
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	50000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	64000
Depth of case	398 mm
Width of case	298 mm
Height of case	215 mm
Case volume	0.0255 cbm
Weight of product	0.0039 kg
Weight of case	4.31 kg
EAN of inner box	4038917105392
EAN case	4038917010542

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-11-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**