



Product Specification
Sample tube, Lithium heparin LH, 4.5 ml, cap green, (LxØ): 75 x 13 mm, with paper label

Page 1



Product description

| | |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number | 32.331.100 |
| Product description | Sample tube, nominal volume: 4.5 ml, cap green, preparation: Lithium heparin, push cap, (LxØ): 75 x 13 mm, round base, with paper label, 100 piece(s)/bag |

Product characteristics

| | |
|----------------------------|------------------|
| Type of preparation | Lithium heparin |
| Preparation concentration | 16 IU/ml blood |
| Colour code | ISO |
| Label/ Print | with paper label |
| Colour of print/label | green |
| Graduation | fill mark |
| Closure type | push cap |
| Base shape | round base |
| Recommended centrifugation | 2000 x g |
| Centrifugation time | 10 Min |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-11-13

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Sample tube, Lithium heparin LH, 4.5 ml, cap green, (LxØ): 75 x 13 mm, with paper label

Page 2

Size

| | |
|----------------------|--------|
| Sample volume | 4.5 ml |
| Diameter | 13 mm |
| Length of product | 75 mm |
| Length excluding cap | 75 mm |

Material & colours

| | |
|-------------------|----------------------------------|
| Colour of product | transparent |
| Closure material | Low Density Polyethylene (LD-PE) |
| Colour of cap | green |
| Tube material | Polypropylene (PP) |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 50000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 48000 |
| Depth of case | 388 mm |
| Width of case | 378 mm |
| Height of case | 215 mm |
| Case volume | 0.0315 cbm |
| Weight of product | 0.0039 kg |
| Weight of case | 4.37 kg |
| EAN of inner box | 4038917105439 |
| EAN case | 4038917010573 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-11-13

This document was prepared by EDP support and is valid without signature.



SARSTEDT