



Product Specification
**Sample tube, Serum CAT, 4 ml, cap black, (LxØ):
75 x 12 mm**

Page 1



Product description

Order number	36.350
Product description	Sample tube, nominal volume: 4 ml, cap black, preparation: Clotting activator, push cap, (LxØ): 75 x 12 mm, round base, 100 piece(s)/bag

Product characteristics

Type of preparation	Clotting activator
Colour code	EU
Graduation	no
Closure type	push cap
Base shape	round base
Recommended centrifugation	2000 x g
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Sample tube, Serum CAT, 4 ml, cap black, (LxØ):
75 x 12 mm

Page 2

Size

Sample volume	4 ml
Diameter	12 mm
Length of product	75 mm
Length excluding cap	75 mm

Material & colours

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	black
Tube material	Polypropylene (PP)

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	2000
Piece(s) / pallet	60000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.0034 kg
Weight of case	7.52 kg
EAN of inner box	4038917081290
EAN case	4038917000437

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT