



### Product description

Order number	41.1506.005
Product description	Micro sample tube Citrate 9NC 0.106 mol/l 3.2%, venous blood collection, preparation: Citrate 9NC, nominal volume: 1.3 ml, (LxØ) with cap: 44 x 10.8 mm, push cap, cap: green, colour code EU, with plastic label, label/print: green, skirted conical base, 100 piece(s)/polystyrene container

### Product characteristics

Type of collection	venous
Type of preparation	Citrate 9NC
Preparation concentration	0.106 mol/l
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	green
Graduation	yes
Cap	flat cap
Closure type	push cap
Application	venous blood collection
Base shape	skirted conical base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Micro sample tube Citrate 9NC 0.106 mol/l 3.2%,  
1.3 ml, push cap, EU**

Volume of work 1,3 ml  
Recommended centrifugation 1500 x g  
Centrifugation time 10 Min

**Size**

Sample volume 1.3 ml  
Diameter 10.8 mm  
Length of product 44 mm  
Length including cap 44 mm  
Length excluding cap 40 mm

**Material & colours**

Colour of product transparent  
Closure material Low Density Polyethylene (LD-PE)  
Colour of cap green  
Tube material Polypropylene (PP)

**Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

**Packaging**

Minimum order qty. 1000  
Type of smallest polystyrene container  
subpackaging  
Piece(s) / inner box 100  
Piece (s) / outer case 1000  
Piece(s) / pallet 44000  
Depth of case 568 mm  
Width of case 348 mm  
Height of case 155 mm  
Case volume 0.0306 cbm  
Weight of product 0.0017 kg  
Weight of case 2.47 kg  
EAN of inner box 4038917129626  
EAN case 4038917022316

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

