



Product description

Order number	01.1602
Product description	S-Monovette® Serum Gel CAT, preparation: Clotting activator/gel, 7.5 ml, membrane screw cap, cap brown, colour code EU/ISO, (LxØ) without cap: 92 x 15 mm, with plastic label, label/print: transparent/brown, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Clotting activator/gel
Colour code	EU/ISO
Label/ Print	with plastic label
Colour of print/label	transparent/brown
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Serum Gel CAT, 7.5 ml, cap brown,
(LxØ): 92 x 15 mm, with plastic label

Page 3

Width of box	134 mm
Height of box	82 mm
Depth of case	440 mm
Width of case	400 mm
Height of case	144 mm
Case volume	0.0253 cbm
Weight of product	0.0088 kg
Weight of case	5.22 kg
EAN of inner box	4038917087209
EAN case	4038917005678

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT