



## Product Specification

**S-Monovette® Lithium heparin gel LH, 7.5 ml, cap orange, (LxØ): 92 x 15 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix**

Page 1



### Product description

Order number	01.1634.088
Product description	S-Monovette® Lithium heparin gel LH, preparation: Lithium heparin gel, 7.5 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 92 x 15 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix 036, label/print: white/orange, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin gel
Preparation concentration	16 IU/ml blood
Colour code	EU
Label/ Print	with plastic label pre-barcoded
Colour of print/label	white/orange
Graduation	fill mark
Cap	membrane screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Lithium heparin gel LH, 7.5 ml, cap orange, (LxØ): 92 x 15 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix**

Page 2

Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	15 Min
Liquid-tight	yes
Connection	membrane screw cap

### Size

Sample volume	7.5 ml
Outer diameter of cap	18 mm
Diameter	15 mm
Length including cap	108 mm
Length excluding cap	92 mm
Length including cap and plunger	129 mm

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	orange
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Lithium heparin gel LH, 7.5 ml, cap orange, (LxØ): 92 x 15 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix**

**Page 3**

### Packaging

Minimum order qty.	20000
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	24000
Depth of box	189 mm
Width of box	134 mm
Height of box	82 mm
Depth of case	440 mm
Width of case	400 mm
Height of case	144 mm
Case volume	0.0253 cbm
Weight of product	0.0088 kg
Weight of case	5.23 kg
EAN of inner box	4038917949682
EAN case	4038917949699

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**