



Product Specification

S-Monovette® Ammonium heparin AH, 9 ml, cap blue, (LxØ): 92 x 16 mm, with paper label

Page 1



Product description

Order number	02.1064.001
Product description	S-Monovette® Ammonium heparin AH, preparation: Ammonium heparin, 9 ml, membrane screw cap, cap blue, colour code EU, (LxØ) without cap: 92 x 16 mm, with paper label, label/print: white, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Ammonium heparin
Preparation concentration	16 IU/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	white
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-11-13

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Ammonium heparin AH, 9 ml, cap blue, (LxØ): 92 x 16 mm, with paper label

Page 3

Width of box	134 mm
Height of box	88 mm
Depth of case	648 mm
Width of case	303 mm
Height of case	150 mm
Case volume	0.0295 cbm
Weight of product	0.0088 kg
Weight of case	5.35 kg
EAN of inner box	4038917125116
EAN case	4038917019767

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-11-13

This document was prepared by EDP support and is valid without signature.



SARSTEDT