



**Product Specification**  
**S-Monovette® EDTA Gel K2E, 9 ml, cap red,**  
**(LxØ): 92 x 16 mm, with paper label**

**Page 1**



### Product description

Order number	02.1333.001
Product description	S-Monovette® EDTA Gel K2E, K2E, preparation: K2 EDTA Gel, 9 ml, membrane screw cap, cap red, colour code EU, (LxØ) without cap: 92 x 16 mm, with paper label, label/print: red, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	K2 EDTA Gel
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	red
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**





**Product Specification**  
**S-Monovette® EDTA Gel K2E, 9 ml, cap red,**  
**(LxØ): 92 x 16 mm, with paper label**

Width of box	134 mm
Height of box	88 mm
Depth of case	648 mm
Width of case	303 mm
Height of case	150 mm
Case volume	0.0295 cbm
Weight of product	0.0093 kg
Weight of case	5.59 kg
EAN of inner box	4038917127677
EAN case	4038917021357

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

